

Renelif Neck Advanced Device



Instructions



Renelif Neck Advanced Skin Therapy

1. Fully charge your RENELIF by connecting it to a USB socket. When fully charged, your RENELIF is ready for use.
2. Turn it on with a long press on the power button. The RENELIF will turn on and a beep will be heard, indicating that normal mode has started.
3. To switch to “cool” mode, press the power button once.
4. To switch to “hot” mode, press the power button again. Turn it off with a long press on the power button.

Your RENELIF has 3 Different Usage Modes:

Normal Mode: Sonic massage therapy helps improve skin vitality, blood circulation, and revitalizes the skin while promoting skin absorption.

Blue LED Light Therapy with Cooling:

The RENELIF will cool down up to -8°C (18°F) below room temperature to help tighten, lift, and treat the skin's appearance.

The blue LED light therapy promotes a healthy skin appearance, reduces the appearance of skin imperfections and impurities, and provides clearer-looking skin. Sonic massage therapy helps improve skin vitality, blood circulation,