

Renelif Eye Skin Therapy

Instructions



Renelif Eye Skin Therapy

Make sure your Renelif eye skin therapy wand is charged before use.

STEP 1. PREPARE YOUR SKIN – Wash and dry your face.

STEP 2. TURN ON YOUR RENELIF EYE SKIN THERAPY WAND – Remove your Renelif eye skin therapy wand from its charging cradle and press the power button to turn it on. You will feel and hear a vibration sound and your Renelif eye skin therapy wand will illuminate blue. You may press the power button again and your Renelif eye skin therapy wand will illuminate blue and begin to warm up. Choose your preferred treatment End.

STEP 3. TREAT YOUR SKIN – Gently glide your Renelif eye skin therapy wand with the metal treatment surface along your under-eye area and your face. The treatment light will automatically turn on upon contact with your skin. Blue light is recommended for treating skin imperfections and should be performed for 3 minutes on each area. Red light is recommended for treating wrinkles and helping with skincare penetration. You may use any skin care with your Renelif eye skin therapy wand and revitalizes the skin while promoting skin absorption.

Red LED Light Therapy with Heating: The RENELIF will heat up to 40°C (104°F). The red LED light therapy treats the appearance