

Voupre Thermal Collection

Instructions



0.25 Derma Roller

Use only after thoroughly cleansing your face. Glide roller across targeted areas, taking care not to apply too much pressure. Continue to roll over the targeted area for 60 seconds at a time. Frequency of use: Every 2 days.

Detoxifying Thermal Mask

On clean, dry skin, apply evenly on face and neck. Leave on for 15 minutes, allowing the warming sensation to fully take effect. Rinse well with water and pat dry.

ENEO Total Device

Use your ENEO Total Device over cream covered areas for 5 minutes. Once complete, remove cream with a warm towel. Once fully removed, wash your face with cold water and apply a moisturizer of choice.

Perfecting Facial Serum

Apply on top of the Detoxifying Thermal Mask, which helps to intensify the mask's effects. To use serum alone, apply all over face and neck, patting into skin until fully absorbed.

Corrective Facial Cream

Apply an even layer on face and neck, massaging this in until fully absorbed.