

Instructions



Pro Collagen Cleavage Décolletage Balm

- We recommend you to use this balm once to twice per week with your Avologi device.
- The balm should be applied as a very small amount onto cleavage and décolletage areas.
- The balm can be applied to your skin either before or after using your Avologi device. Both ways are extremely effective. Make sure to leave the balm on your skin once you've completed your session.



Thermikos Stretch Marks Age Defying Concentrate

- We recommend you use this treatment at least once per month in areas below the neck.
- Mix the mask before each use.
- Apply a small amount of the concentrate to desired areas and incorporate a 5 minute session with any of your Avologi devices.
- After 5 minutes of treatment, remove the concentrate with a warm towel.
- Once complete, make sure to tighten the jar and store until next to use.



OTC Age Defying Treatment Gel

- Apply a smooth layer of the gel on to the treatment area of your body before using your Avologi device.
- Turning on your device, glide it gently in a slow upward motion.
- For optimal results, use your Avologi device for five minutes per area, and at least three times per week for the first Month. Following this you may use it twice per month for maintenance.
- Once complete, massage excess gel into skin and clean your device with a dry cloth. Suitable for all skin types and tones.